

Zero-Gravity vs. Conventional Lead Shields

Comparison of operator eye exposures when working from femoral region, side or head of patient.⁷

HEAD

215 Fluoroscopy
Minutes

0.022 Operator Exposures
(microSV/min)

SIDE


67 Fluoroscopy
Minutes

0.035 Operator Exposures
(microSV/min)

FEMORAL

47 Fluoroscopy
Minutes

0.013 Operator Exposures
(microSV/min)



HEAD

97 Fluoroscopy
Minutes

1.817 Operator Exposures
(microSV/min)

SIDE

32 Fluoroscopy
Minutes

4.762 Operator Exposures
(microSV/min)

FEMORAL

35 Fluoroscopy
Minutes

3.466 Operator Exposures
(microSV/min)