



Navi-Crani™

Maintain sterility during Navigated Craniotomy procedures



Introducing Navi-Crani™ - the sterile patient drape designed to reduce the risk of contamination during navigated craniotomy procedures

AORN Guidelines state:

“Perforation of barrier materials may provide portals of entry and exit for microorganisms, blood and other potentially infectious materials.”

Navi-Crani is installed without the need to cut a hole in the drape to fit over the navigation system arm, thereby decreasing potential contamination.

Navi-Crani works well with Brainlab, Stryker and Medtronic navigation systems.

Permitting a quicker install, Navi-Crani leads to a faster room turnaround



PATENT PENDING

Visit CFImedical.com to find out more information.

Part #5568
Box of 8

Call 810-750-5300 to learn more about purchase options.

*Perioperative Standards and Recommended Practices, 2014 Edition; Publisher: AORN, Inc.; ISBN: 978-1-888460-85-8

NAVI-CRANI™

a CFI MEDICAL System

70004-08-002

IMPROVING EFFICACY